Scientific Abstracts

12 years of clinical research | **Vol. 2**

Dear Doctor,

A rolling stone gathers no moss.

What is valid for stones, and certainly also for human beings, is probably also a good guideline for scientific research. This is why two years after vol. 1 of mectron PIEZOSURGERY®'s abstract book, we are now pleased to publish vol. 2.

The advantages of piezoelectric bone surgery have led to a continuously higher acceptance in daily clinical practice and, of course, more research. Therefore, in addition to various chapters on oral and dental surgery, in vol. 2 you will also find a few abstracts of publications on maxillofacial applications. This will give an idea of the extraordinary potential of PIEZOSURGERY® technology.

Evidence based knowledge remains the foundation from which all current and future steps in the evolution of mectron PIEZOSURGERY® are made. We hope that you as an expert in bone surgery will appreciate the present collection as a proof statement for this philosophy.

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Muñoz-Guerra MF, Naval-Gías L, Capote-Moreno A.

Le Fort I osteotomy, bilateral sinus lift, and inlay bone-grafting for reconstruction in the severely atrophic maxilla: a new vision of the sandwich technique, using bone scrapers and PIEZOSURGERY®.

ENGLISH J Oral Maxillofac Surg. 2009 Mar;67(3):613-8.

■ **Le Fort I osteotomy, bilateral sinus lift, and inlay bone-grafting for reconstruction in the severely atrophic maxilla: a new vision of the sandwich technique, using bone scrapers and PIEZOSURGERY®.**

Severe atrophy of the edentulous maxilla and progressive pneumatization of the maxillary sinus can compromise the insertion of dental implants. In this context, ideal implant positioning is limited by inadequate height, width, and quality of the bone. Le Fort I osteotomy and interpositional bone graft is an excellent treatment concept for the dental rehabilitation of patients with atrophied maxilla and reversed intermaxillary relationship. In this report, we indicate the transcendent aspect of elevation and preservation of maxillary sinus and nasal mucosa, modifying the sandwich technique by the useful of bone scrapers and PIEZOSURGERY®. The procedure is described including a 1-stage approach using cortico-cancellous bone blocks through which implants are placed. In the extremely atrophied alveolar process of the maxilla, this technique provides the desired gain of bone, allows for the ideal placement of dental implants, and improves any discrepancy between the upper and lower arches.

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Palatal versus vestibular piezoelectric window osteotomy for maxillary sinus elevation: a comparative clinical study of two surgical techniques.

ENGLISH Oral Surg Oral Med Oral Pathol Oral Radiol Endod. 2009 May;107(5):648-55.

■ **Palatal versus vestibular piezoelectric window osteotomy for maxillary sinus elevation: a comparative clinical study of two surgical techniques.**

Objectives: The goal of this study was to compare the surgical advantages and disadvantages of a new palatal access osteotomy for sinus elevation with a conventional lateral approach.

Study design: In 32 patients, either a palatal (n = 16) or a lateral (n = 16) osteotomy to the maxillary sinus was performed under local anesthesia. The palatal access included a circular paramarginal incision and elevation of a palatal mucosal flap based on a median pedicle. The lateral access was performed by

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vestibular standard incision and development of a mucoperiosteal flap with a vestibular and superior basis. For all osteotomies a piezoelectric device was used. The sinus cavity was augmented with synthetic nanostructured hydroxyapatite graft material.

Results: Intraoperative complications during both procedures were minimal and wound healing was uneventful. Membrane perforation occurred in 19% of the palatal group and in 19% of the lateral group. Soft tissue management of the palatal technique was superior to that of the lateral approach, because the vestibular anatomy was not altered and consequently no disharmonious soft tissue scarring and no postoperative swelling occurred.

Conclusion: The palatal approach permitted higher postoperative comfort, especially for edentulous patients, because full dentures could be incorporated directly after surgery with almost perfect fit.

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■ Augmentation of alveolar ridge defects: autologous bone transplant from the zygomatic alveolar crest--a new technique

The insertion of dental implants with simultaneous loss of alveolar bone still represents a challenge. A traumatic tooth loss within the range of the front of the maxillary is frequently connected with a clear substance defect. The goal of a pre-implantologically treatment must be therefore the hard tissue reconstruction of the defect. Thus the implant

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camp is optimized and the prosthetical necessary implant position is guaranteed. In this case report the reconstruction of a vertical and horizontal bone loss with a transplant from the Crista zygomatico alveolaris is represented. The withdrawal of the bone transplant took place under indulgence of the vestibular membrane of the maxillary sinus by means of PIEZOSURGERY®. After an appropriate healing time the single tooth gap was supplied with an ITI implant of the company Straumann (Freiburg, Germany).

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Gellrich NC, Held U, Schoen R, Pailing T, Schramm A, Bormann KH.

Alveolar zygomatic buttress: A new donor site for limited preimplant augmentation procedures.

ENGLISH

J Oral Maxillofac Surg. 2007 Feb;65(2):275-80.

■ Alveolar zygomatic buttress: A new donor site for limited preimplant augmentation procedures.

Purpose: Augmentation of anterior maxillary alveolar bone defects before placement of endosseous implants still poses a clinical challenge. Previous traumatic loss of upper anterior teeth is often concomitant with significant bone loss. The aim of augmentation is to reconstitute as far as possible the original hard tissue contour. This is a prerequisite for later optimal positioning of the dental implants in accordance with prosthodontic and functional principles.

Materials and Methods: This article describes a new method for the reconstruction of small anterior maxillary alveolar bone defects using donor bone from the zygomatic buttress region.

Results: This region provides harvested bone with a natural convex shape, which is ideally suited for the anterior alveolar process region. Later fibrous tissue transplantation to provide more anterior projection to the anterior dental area is also avoided.

Conclusion: Bone harvesting was done with the piezoelectric device, which avoids trauma to the mucosal membrane of the maxillary sinus.

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■ Piezo surgery and bone grafting

Piezo surgery significantly simplifies graft harvesting from the ramus, allowing a bone graft to be taken in a site where the postoperative course is always straightforward without a risk of neurological sequela and with mucosal healing of a quality not obtained with harvesting from the symphysis.

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Relocation of a malpositioned maxillary implant with piezoelectric osteotomies: a case report.

ENGLISH

Int J Periodontics Restorative Dent. 2008 Oct;28(5):489-95.

■ Relocation of a malpositioned maxillary implant with piezoelectric osteotomies: a case report.

Implant relocation is a new surgical technique for correcting the alignment of malpositioned implants by mobilizing them with the surrounding bone until the desired position is achieved. In this case report, a 25-year-old woman was treated for the malposition of an implant in the maxillary left canine site. The use of a piezoelectric scalpel permits narrow, precise, and safe osteotomies, thus preventing involvement of the soft tissue and producing better healing potential compared to burs or saws. The results suggest that inadequately axially inclined implants can be successfully reconfigured using segmental piezoelectric osteotomies.

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■ **Juxta-canal implants**

Implantology forms part of practitioners' therapeutic arsenal and can be integrated in treatment plans for oral rehabilitation. It obviates the need for removable or fixed prostheses, which sometimes necessitate preparation of healthy teeth. Nevertheless, implant placement depends on the presence of sufficient bone volume at the implant site. If necessary, prior performance of a bone graft may have to be planned (1). Another solution is maximum use of the available bone. The presented clinical case describes an operative technique that allows use of the entire bone height available above the inferior alveolar nerve in the posterior area.

This technique employs piezoelectric bone surgery, developed by Vercellotti (11), an interesting alternative to the classical techniques of bone surgery (2). Piezo surgery has been conceived to simplify surgical protocols and improve the predictability of the results and the safety of the surgical procedure.

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ENGLISH Oral Surg Oral Med Oral Pathol Oral Radiol Endod. 2009 Dec;108(6):e1-5.

■ **Bone lid technique with PIEZOSURGERY® to preserve inferior alveolar nerve.**

The inferior alveolar nerve is at risk during surgical removal of impacted mandibular teeth and excision of benign tumors of mandible. Manual and/or mechanical instruments used in the close proximity of delicate structures do not allow the control of the cutting depth and can damage these structures by accidental contact. Piezoelectric surgery is a new and innovative bone surgery technique using the microvibrations of special scalpels at ultrasonic frequency so therefore soft tissue will not be

■ **Tecnica del „bone lid“ con PIEZOSURGERY® per la preservazione del nervo alveolare inferiore.**

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damaged even upon accidental contact with the cutting tip. This article presents an alternative technique that uses piezoelectricity to minimize trauma to the inferior alveolar nerve, vascular tissues, or surrounding dental tissues.

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■■■ EXTRACTIONS

Sortino F, Pedullà E, Masoli V.

The piezoelectric and rotatory osteotomy technique in impacted third molar surgery: comparison of postoperative recovery.

ENGLISH

J Oral Maxillofac Surg. 2008 Dec;66(12):2444-8.

■ **The piezoelectric and rotatory osteotomy technique in impacted third molar surgery: comparison of postoperative recovery.**

Purpose: The aim of this study was the comparison of postoperative outcome in mandibular impacted third molars treated by piezoelectric surgery or by rotatory osteotomy technique.

Patients and methods: One hundred patients with impacted mandibular third molars were included in the study. Fifty patients were treated by rotatory osteotomy technique (group A) and 50 patients were treated by piezoelectric osteotomy technique (group B). Therapeutic protocol was the same for both groups. Twenty-four hours after surgery, 2 different parameters, facial swelling and trismus, were evaluated in both groups. A pair of compasses were used for the evaluation of facial swelling and trismus was evaluated.

Results: The average surgery time was 17 minutes in group A, 23 minutes in group B; the mean facial swelling was 7.04 mm in group A, 4.22 mm in group B; trismus was 16.76 mm in group A, 12.52 mm in group B. Statistical analysis showed a significant reduction ($P < .05$) of postoperative facial swelling and trismus in group B; however, in this group, a statistically significant increased ($P < .05$ vs group A) surgery time was required.

Conclusion: The piezoelectric osteotomy technique produced a reduced amount of facial swelling and trismus 24 hours after surgery, but a longer surgery time was required when compared with the rotatory osteotomy technique.

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Seshan H, Konuganti K, Zope S.

PIEZOSURGERY® in periodontology and oral implantology

ENGLISH

J Indian Soc Periodontol. 2009 Sep;13(3):155-6.

■ PIEZOSURGERY® in periodontology and oral implantology

Abstract

Periodontitis is a chronic inflammatory disease of the tooth-supporting structures. The treatment of this condition is largely based on the removal of local factors and restoration of the bony architecture. Moreover, in the era of modern dentistry, successful implant therapy often requires sound osseous support. Traditionally, osseous surgery has been performed by either manual or motor-driven instruments. However, both these methods have their own advantages and disadvantages. Recently, a novel surgical approach using piezoelectric device has been introduced in the field of periodontology and oral implantology. This article discusses about the wide range of application of this novel technique in periodontology.

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The reconstruction of a maxillary anterior dento-alveolar defect in patients with trauma has been a challenge for surgeons. Extensive loss of bone and teeth in the anterior maxilla presents a complex problem for reconstruction. This is owing to the difficulty in achieving complete closure using overlying soft tissue. Tension-free sutures cannot be used after a large bone graft because the overlying soft tissue on severe bone defects of the anterior maxilla is often deficient and is attached to the underlying atrophic bone by scarring. Distraction osteogenesis provides a method to regain both hard tissue and soft tissue without any grafting. We describe a patient who had severe maxillary anterior bony defects that were restored by means of piezoelectric distraction osteogenesis, followed by dental implant placement. Clinical, radiological, and histological results showed that the reconstruction was successful.

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Stübinger S, Goethe JW.

Bone Healing After PIEZOSURGERY® and its Influence on Clinical Applications

ENGLISH

Journal of Oral and Maxillofacial Surgery 2007 Sep;65(9):39.e7-39.e8.

■ Bone Healing After PIEZOSURGERY® and its Influence on Clinical Applications

Statement: Osteotomies are well-established procedures in oral and maxillofacial surgery for either corrective measures or for gaining access to subjacent hard and soft tissues. Yet surgical reconstruction of bony defects can often be very challenging as thin and fragile bony structures are especially prone to fracture due to bulky cutting tips or massive application of pressure by conventional mechanical instruments. By employing the new technique of PIEZOSURGERY® such complications can be overcome. It was the aim of this animal study to analyse the bone remodelling process after ultrasound osteotomy in the sheep tibia under functional loading after 8 and 12 weeks and to evaluate its impact on clinical applications.

Materials and Methods: In 12 sheep the PIEZOSURGERY® device (Mectron, Italy) was used to perform a midshaft osteotomy of the tibia diaphysis. The instrument operates with modulated ultrasound and thereby generates micromovements between 60 and 200 mm/sec. Physiological sodium chloride solution of approx. 4°C was used as an irrigant. For all osteotomies bone saw OT7 with parameters mode boosted burst c and pump 5 were used according to the manufacturer's recommendation. All osteotomies were performed by one surgeon in a latero-lateral movement as well as in a turning guidance around the bone.

Method of Data Analysis: After 2 and 3 months specimen were evaluated by radiographic and histologic analysis (toluidin-blue-staining, fuchsin-vitalstaining). A complete, semi-quantitative and histomorphometrical evaluation was undertaken including fluorochromatic stainings indicative for bone remodeling. Data were compared to results gained by analogous experiments performed with other osteotomy techniques.

Results: Piezoelectric osteotomy permitted a micrometric, selective cut and a clear surgical site by the cavitation effect created by irrigation/cooling solution and oscillating tip. No excessive bleeding was encountered. The evaluation of the PIEZOSURGERY® specimen proved 8 weeks as well as 12 weeks after surgery an ingrowth of vital bone-forming tissue into the osteotomy gap. The remodeling in the compact bone was undisturbed and the osteotomy gap was filled with new bone. Additional radiological findings supported these findings. The bone fragments were completely healed, the bone marrow cavity restored as well as the external callus formation was subsided undergoing piezoelectric

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surgery. In general bone healing was faster than known from conventional methods.

Conclusion: PIEZOSURGERY® definitely enhances handling of delicate structures in the oral and maxillofacial region. Concerning osteotomies of thin and fragile bones, application of ultrasound is assessed to be superior to other mechanical instruments, because of easy handling, efficient bone ablation and minimal accidental harm to adjacent soft tissue structures. As bone healing is not disturbed by the PIEZOSURGERY®, but even seems to be improved, this method will have a major influence on new minimally invasive bone surgery techniques with special regard to biomechanics.

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Sohn DS, Lee JK, An KM, Shin HI.

Histomorphometric evaluation of mineralized cancellous allograft in the maxillary sinus augmentation: a 4 case report.

ENGLISH

Implant Dent. 2009 Apr;18(2):172-81.

■ Histomorphometric evaluation of mineralized cancellous allograft in the maxillary sinus augmentation: a 4 case report.

Purpose: The aim of this article is to report the histomorphometric analysis and clinical study of mineralized cancellous bone allografts in maxillary sinus augmentation procedures in 4 cases.

Materials and Methods: Sinus bone augmentation, using a piezoelectric device and mineralized cancellous allografts, was performed in 3 patients. After an average of about 8 months of healing, a trephine core sample was gained from the lateral windows area for histomorphometric analysis. The mineralized cancellous allografts demonstrated normal healing of hard tissues.

Results: Upon reentry for dental implant placement or secondary surgery, the graft materials showed hard bony structures within the maxillary sinus. Histomorphometric analysis of the cores revealed an average vital new bone content of 17.88%.

Conclusion: Mineralized cancellous allografts showed favorable bone regeneration histologically and was clinically useful in augmenting bone volume for implant placement in the sinus bone graft.

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■ Sound levels of the PIEZOSURGERY®. Risk of permanent damage to hearing

In the past, research has regularly been carried out concerning the sound levels of various drilling devices and the impact these have on those who regularly use these devices. The present research is concerned with the possible permanent damage to hearing which can occur during the use of a newly developed instrument, the PIEZOSURGERY®. Measurements have been performed in a standardized set-up in which use is made of the lower jaws of pigs and a weight scale for measuring various degrees of pressure. The boundary values of the permissible exposure to noise were determined. The values of the PIEZOSURGERY® were compared with conventional drilling machines. It was concluded that using the PIEZOSURGERY® for less than 1.5 hours per day implies no risk of permanent damage to hearing. This means that the risk in the case of PIEZOSURGERY® is somewhat lower than that of conventional drilling devices.

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■ **Light microscopic examination of rabbit skulls following conventional and ultrasonic osteotomy**

Introduction: The novel ultrasonic osteotomy technique (PIEZOSURGERY®) is an alternative to conventional osteotomy devices. The aim of the present study was to carry out morphological comparison of the bone surface using conventional osteotomy techniques in comparison to the rather new ultrasonic osteotomy technique by means of a reflected-light microscopic examination.

Materials and Methods: Following the sacrifice of 12 rabbits, 24 standardized bone samples were removed from the skull. The osteotomy devices used were a rotating instrument (Lindemann bur), an oscillating micro-saw, and an ultrasonic osteotomy device (PIEZOSURGERY®) with insert tips OT6 and OT7. The times needed for osteotomy were measured. The bone surfaces were examined using a reflected-light microscope with a magnification of 40x and 100x.

Results: Osteotomy with PIEZOSURGERY® is significantly more time consuming than osteotomy with conventional methods ($p < 0.05$). Following osteotomy with the ultrasonic device, the reflected-light microscopic examinations of the unmodified bone samples revealed typical bone structure of the calvaria, including compacta externa, diploe and compacta interna. On the contrary, following osteotomy with the conventional devices, the diploe structure presented distinct modifications. The cancellous spaces were filled with bone debris, and the cancellous structure was demolished. The samples prepared by the micro-saw technique showed a superficially condensed and grooved surface.

Conclusion: In the present study, well-defined differences were observed following osteotomy with conventional devices and osteotomy with the ultrasonic device. The integrity of the bony structure observed after the ultrasonic technique could benefit the bone healing process. Further studies dealing with the bone healing process after using different osteotomy techniques are recommended.

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ITALIANO

Minerva Stomatol. 2006 Apr;55(4):199-207.

■ Peri-implant tissue healing in implant placed in post-extraction sockets of maxillary premolars.

Aim: The aim of this study was the clinical and radiographic evaluation of peri-implant tissues healing associated with two-stage implants performed with one-stage surgery and placed into post-extraction sockets of maxillary premolars.

Methods: Ten natural tapered submerged titanium implants with double acid-etched surface were installed in post-extraction sockets of maxillary premolars showing no bone dehiscence or fenestrations. Implants were selected according to the greatest dimension compatible to vertical and horizontal diameters of the post-extraction sockets. Peri-implant defects showing more than 2 mm of distance between the marginal palatal bone and the implant margin, were not treated with a regenerative procedure but with a bone swaging technique by means of piezoelectric surgery. At 16 weeks osteointegration of implants was examined with a clinical and standardized radiographic evaluation.

Results: All patients showed good clinical healing without any complication. At 16 weeks all implants, either those with reduced marginal defect or those with modified crestal bone profile, were osteointegrated.

Conclusion: The use of two-stage implants performed with one-stage surgery in 4 walls post-extraction sockets of maxillary premolars, seems to allow osteointegration both in reduced peri-implant horizontal defects as well as in defects larger than 2 mm. Furthermore, these procedures seem to simplify the management of soft peri-implant tissues.

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■ The ultrasonic cut in implantology; 1st part: technologies

The ultrasonic lancet is a surgical device able to cut out with precision hard tissue and to facilitate the solid interfaces cleavage. It uses microvibrations of intermediate frequency generated by a piezoelectric transducer and applied to titanium nitride hardened and/or diamond-coated inserts. With its vast range of inserts, it finds many applications in oral and maxillofacial surgery, such as non-traumatic dental avulsions, root planing and bone defect debridement, or cyst removal. It also proposes a simplified protocol for the sinus lift surgery. Lastly, it offers a true revolution in the bone grafting surgery by allowing precise and non-traumatic graft harvesting. This article presents its physical, technological and clinical aspects and discusses its most promising applications. Indeed, although the ultrasonic lancet remains a safe and powerful tool in many circumstances, it is however significant to define a reasoned application field to it.

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■ **The ultrasonic cut in implantology; 2nd part: clinical applications**

The ultrasonic lancet makes possible to realize with extreme precision and safety a great number of oral and maxillofacial surgeries. However, it is during four quite particular interventions that this tool will present really higher performances than the conventional instrumentation: non-traumatic removal of osteo-integrated implants, symphysis and retromolar bone harvesting and inferior alveolar nerve lateralization. The objective of this article is to present in detail each one of these four applications and to discuss their advantages and disadvantages compared to the former techniques. Indeed, the use of the piezoelectric lancet considerably simplifies these surgical protocols with non-traumatic cut and ultrasonic solid interfaces cleavage, which supports the bone fragments rupture.

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Within a few years of ultrasonic technology, a very powerful ultrasonic device has been developed. This device, together with its modified tips, is a tool that allows impressive cutting, providing numerous very practical solutions in implant surgery. This review article describes the principles of piezoelectric osteotomy as well as a survey of the numerous known clinical applications of this instrument in implant surgery.

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■ The piezoelectric scalpel: a delicate tool highly suitable for use in dentistry

Bone surgery (osteotomy and osteoplasty) has gained in refinement through PIEZOSURGERY®. The principle is based on use of an ultrasonic scalpel with inserts connected to a handpiece identical with that of a scaler, allowing a very fine cut, up to 0.2 mm for the finest inserts.

The effect on mineralized tissues is obtained by means of double oscillation of the inserts, a vertical amplitude of 20 to 60 µm and a horizontal amplitude of 60 to 200 µm, by using ultrasonic frequencies that can be modulated from 22 to 30 kHz and simultaneously combining these frequencies.

In this frequency range, the second advantage is the absence of any action of the inserts on the soft tissues, which is of particular interest with regard to anatomical structures at risk (neurovascular bundles, mucosa).

We present PIEZOSURGERY® by means of five clinical cases showing some of the possibilities of use of the system.

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Robiony M, Polini F, Costa F, Vercellotti T, Politi M.

Piezoelectric Bone Cutting in multipiece maxillary osteotomies.

Technical Note.

ENGLISH

J Oral Maxillofac Surg. 2004; 62:759-761.

■ Piezoelectric Bone Cutting in multipiece maxillary osteotomies. Technical Note.

PIEZOSURGERY® (patented by Mectron Medical Technology, Carasco [Genova], Italy.) is a new innovating technique used to perform safe and effective osteotomies using piezoelectric ultrasonic vibrations. It was first reported for preprosthetic surgery, alveolar crest expansion, and sinus grafting. We introduce and report the use of PIEZOSURGERY® for multipiece maxillary osteotomies, to overcome many of the complications of this delicate surgery on hard and soft tissues. Because of its micrometric and selective cut, the piezoelectric device produces safe and precise osteotomies without any osteonecrosis damage. This device works only on mineralized tissues, sparing soft tissues and their blood supply. Among various surgical phases, osteotomy is one of the most technique-sensitive procedures in maxillofacial surgery. Osteotomies are usually conducted close to delicate anatomic structures, such as vestibular and lingual/palatal soft tissues that provide bone vascularization through the periosteum. Furthermore, bone is a hard tissue and many cutting or drilling osteotomes are very crude tools. In particular, rotating instruments are potentially injurious, due to the production of excessively high temperatures during osseous drilling, which can produce marginal osteonecrosis and impair bony regeneration. It is widely and strongly recommended to use a careful surgical technique and to reduce the amount of frictional heating with saline solution irrigation. Among cutting techniques, PIEZOSURGERY® is a new and innovating method that uses piezoelectric ultrasonic vibrations to perform precise and safe osteotomies, due to its characteristics of a micrometric and selective cut. It was first invented by Tomaso Vercellotti, MD, DDS, to overcome the limits of traditional instruments in oral bone surgery. The purpose of this report was to present the use of the piezoelectric cut in segmental maxillary Le Fort I osteotomy, a field in which effectiveness, precision, and safety of osteotomies are of paramount importance.

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Eggers G, Klein J, Blank J, Hassfeld S

PIEZOSURGERY®: an ultrasound device for cutting bone and its use and limitations in maxillofacial surgery.

ENGLISH

Br J Oral Maxillofac Surg. 2004 Apr;42(5):451-3.

■ **PIEZOSURGERY®: an ultrasound device for cutting bone and its use and limitations in maxillofacial surgery.**

PIEZOSURGERY® uses modulated ultrasonic vibration to allow controlled cutting of bony structures. Delicate bony structures can be cut easily and with great precision, without destruction of soft tissue. We have found this device useful when exact cutting of thin bones is essential. However, it is of only limited use in cutting thick bones and in regions with limited access.

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■ **What is PIEZOSURGERY®, Two-years experience in craniomaxillofacial surgery**

Introduction: PIEZOSURGERY® is a new surgical technique used in dentistry to section hard tissues without damaging adjacent soft tissues. We hypothesized that such a device could also be useful in craniofacial and orthognathic surgery.

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Material and method: An ultrasonic device (Mectron) was employed in different craniofacial surgical procedures: a) to perform 144 Le Fort I osteotomies, 140 palatal expansions after Le Fort I osteotomies, and 134 bilateral sagittal osteotomies; b) to perform a Le Fort III osteotomy for treatment of Crouzon syndrome in 2 patients; c) to perform 5 segmental osteotomies and 3 osteotomies of the inferior edge of the mandible for facial asymmetry; d) to perform 12 cases of unicortical calvarial bone grafting; to remove the superior orbital roof in 20 cases of craniofaciostenosis and the frontal bone in 5 cases; e) to remove the external wall of the orbit or the anterior and posterior wall of the frontal sinuses in 10 cases of orbital cavity tumors; f) to approach the skull base through the frontal sinuses in 4 cases.

Results: Analysis of the results showed that PIEZOSURGERY®: a) allows very precise cutting; avoids bone cutting using an osteotome; b) spares soft tissue such as brain, dura-mater, palatal mucosa; c) preserves the nerves, especially those of the inferior alveolar; d) increases the time of bone cutting but not the overall operative time because of the absence of soft tissue protection.

Discussion: PIEZOSURGERY® is a new technical procedure, which can be particularly advantageous for bone cutting in maxillo facial surgery, sparing adjacent soft tissues such as brain, palatal mucosa, and the inferior alveolar nerve from any damage. The device's lack of power appears to be a minor problem compared with the advantages.

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Beziat JL, Béra,JC, Lavandier B, Gleizal A.

Ultrasonic osteotomy as a new technique in craniomaxillofacial surgery.

ENGLISH

International Journal of Maxillo-facial Surgery, 2007;36(6):493-500.

■ Ultrasonic osteotomy as a new technique in craniomaxillofacial surgery.

Ultrasound osteotomy is a new surgical technique used in dentistry to section hard tissues without damaging adjacent soft tissues. It was hypothesized that this could also be useful in craniofacial and orthognathic surgery. An ultrasonic device was employed in the following craniofacial surgical procedures: 144 Le Fort I osteotomies, 140 palatal expansions after Le Fort I osteotomies and 140 bilateral sagittal osteotomies; 2 Le Fort III osteotomies for treatment of Crouzon syndrome in two patients; 12 cases of unicortical calvarial bone grafting; removal of superior orbital roof in 25 cases of craniofaciostenosis; removal of external wall of the orbit in 10 cases of orbital cavity tumour; removal of anterior and posterior walls of the frontal sinuses in four cases of orbital cavity tumour. Integrity of soft

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tissues and surgical time were evaluated. Functional results were good without any soft-tissue damage being observed, but the overall operative time was increased. Ultrasound osteotomy is a new technical procedure that is advantageous for bone cutting in multiple situations, with minimal to no damage in adjacent soft tissues such as brain, palatal mucosa and the inferior alveolar nerve.

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Gleizal A, Béra JC, Lavandier B., Béziat JL.

Piezoelectric osteotomy: a new technique for bone surgery – advantages in craniofacial surgery.

ENGLISH

Childs Nerv Syst . 2007;23(5):509-513.

■ Piezoelectric osteotomy: a new technique for bone surgery – advantages in craniofacial surgery.

Introduction: Ultrasonic bone cutting is a new surgical technique used in dentistry to section hard tissues without damaging adjacent soft tissues.

We hypothesized that such a device could also be useful in craniofacial surgery, particularly during the removal of the superior orbital roof during cranio-faciostenosis surgery.

Materials and Methods: An ultrasonic device was employed in different craniofacial surgical procedures: 1. to remove the superior orbital roof in 30 cases of craniofaciostenosis, 2. to perform a Le Fort III osteotomy for the treatment of Crouzon syndrome in two patients, 3. to cut the parietal and frontal bone in 30 cases of craniofaciostenosis. The integrity of soft tissues and surgical time was evaluated.

Results: Functional results were good without any soft tissue damage appreciated. The overall operative time, however, was increased.

Conclusions: PIEZOSURGERY® is a new technical procedure, which can be advantageous for bone cutting in multiple situations with minimal to no damage in adjacent soft tissues.

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Guo ZZ, Liu X, Li Y, Deng YF, Wang Y.

The use of PIEZOSURGERY® osteotomy in treatment of longstanding maxillary fractures: report of 12 consecutive patients.

ENGLISH

Shanghai Kou Qiang Yi Xue. 2007 Feb;16(1):97-9.

■ The use of PIEZOSURGERY® osteotomy in treatment of longstanding maxillary fractures: report of 12 consecutive patients.

Purpose: To evaluate the clinical applicability of PIEZOSURGERY® osteotomy: a new safe technique in managing long standing maxillary fractures.

Methods: 12 patients with long-standing maxillary fractures were surgically treated using Le Fort I osteotomy. During operation, PIEZOSURGERY® osteotomy was used for bone cutting and splitting. After repositioning, the bone segments were rigidly fixed with micro Ti-plate, Ti-mesh. All the patients were followed up for 6 to 12 months, and the functional and esthetic results were evaluated.

Results: Ultrasonic microvibrations allow accurate bone cutting without oscillating injuries to the soft tissue. All the wounds healed primarily without complications. The postoperative occlusion and appearance were satisfactory.

Conclusions: Maximal recovery of mastication and appearance can be achieved by using PIEZOSURGERY® osteotomy with fixation materials such as Ti-plates and Ti-meshes in selected patients with long-standing maxillary fractures.

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ITALIANO

G Chir 2007 May; 28(5):222-6.

■ Mandibular enostosis removal with PIEZOSURGERY®: case report.

The bone surgery has always used manual and rotary instruments. The biomedical engineering with ultrasound working principle has given a new surgery instruments, which allow a selective cutting action of bone tissue and the protection of soft tissue. Our case shows an uncommon clinical lesion surgically dangerous for the narrow adjoining of important anatomical structures as the lower alveolar artery and the inferior alveolar nerve. The clinical result and recovery time go toward a smaller traumatic situation of this methodology of the cutting of bone tissue.

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Robiony M, Polini F, Costa F, Zerman N, Politi M.

Ultrasound bone cutting for surgically assisted rapid maxillary expansion (SARME) under local anaesthesia.

ENGLISH

Int J Oral Maxillofac Surg. 2007;36(3):267-9.

■ Ultrasound bone cutting for surgically assisted rapid maxillary expansion (SARME) under local anaesthesia.

Surgically assisted rapid maxillary expansion (SARME) is a well-established therapy for correction of maxillary transverse deficiency in adults, when consolidation of sutures has just been completed. It can be performed either under general or under local anesthesia and it can be accomplished with many surgical techniques. One of the most critical steps of SARME is the detachment of the pterygo-maxillary junction, due to the risks connected to such procedure. When required to obtain specific expansion patterns, the pterygo-maxillary separation has been suggested until now only for interventions under general anesthesia, due to the dangerousness and the rawness of this surgical step in awake patients. The authors introduce the use of an ultrasonic bone-cutting device to perform all osteotomic steps of SARME under local anesthesia on an outpatient basis, including pterygo-maxillary detachment. This ultrasonic device is unique in that the osteotomic action occurs only when the tool is employed on mineralized tissues, while it stops on soft tissues. It works in a linear pattern of vibration and it allows precise osteotomies without producing any heat damage to osteotomic surfaces and without any dangerous hammer-related stroke. Due to its precision and safety, this device named PIEZOSURGERY®, allows patients to undergo all the steps of SARME under local anesthesia, also without hospitalization.

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Critical evaluation of piezoelectric osteotomy in orthognathic surgery: operative technique, blood loss, time requirement, nerve and vessel integrity.

ENGLISH

J Oral Maxillofac Surg. 2008 Apr;66(4):657-74.

■ Critical evaluation of piezoelectric osteotomy in orthognathic surgery: operative technique, blood loss, time requirement, nerve and vessel integrity.

Purpose: Piezo-osteotomy feasibility as a substitute for the conventional saw in orthognathic surgery was evaluated regarding operative technique, blood loss, time requirement, and nerve and vessel integrity.

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Patients and methods: Fifty patients had orthognathic surgery procedures in typical distribution using piezosurgical osteotomy: 22 (44%) monosegment, 26 (52%) segmented Le Fort I osteotomies; 48 (48%) sagittal split osteotomies, 6 (12%) symphyseal, and 4 (4%) mandibular body osteotomies. Controls were 86 patients with conventional saw and chisel osteotomies: 57 (66%) monosegment, 25 (29%) segmented Le Fort I osteotomies, 126 (73%) sagittal split, and 4 (5%) symphyseal osteotomies.

Results: Piezosurgical bone osteotomy permitted individualized cut designs, enabling segment interdigitation after repositioning. Angulated tools weakened the pterygomaxillary suture; auxiliary chisels were required in 100% of cases for the nasal septum, and lateral nasal walls as 46% pterygoid processes. After downfracture, the dorsal maxillary sinus wall and pterygoid processes were easily reduced. Hemorrhage was successfully avoided with average blood loss of 541 +/- 150 mL versus 773 +/- 344 mL (P = .001) for a conventional bimaxillary procedure. Sagittal mandibular osteotomy required considerable time (auxiliary saw in 13%); the lingual dorsal osteotomy was mostly performed tactile. Time investment remained unchanged: 227 +/- 73 minutes per bimaxillary standard osteotomy versus 238 +/- 61 minutes (P = .5); clinical courses and reossification were unobtrusive. Alveolar inferior nerve sensitivity was retained in 95% of the study collective versus 85% in the controls (P = .0003) at 3 months postoperative testing.

Conclusions: Piezoelectric osteotomy reduced blood loss and inferior alveolar nerve injury at no extra time investment. Single cases require auxiliary chiseling or sawing. Piezoelectric drilling for screw insertion and complex osteotomy designs may be developed to maintain bone contact or interdigitation after repositioning and minimize need for osteofixation.

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■ PIEZOSURGERY® for genioglossal advancement in treatment of obstructive sleep apnea syndrome. Summary

The tongue is a frequent cause of upper airway obstruction in obstructive sleep apnea syndrome. The genioglossi muscles are the front-most muscles of the tongue, attached to the geniotubercle on the inner cortical of the mandibular symphysis. The collapsibility of the tongue can be managed with genioglossus advancement and the obstructive sleep apnea syndrome decreased. Piezo surgery can decrease the risk of apical dental lesion. The perfect depth control during osteotomy and absence of soft tissue lesion considerably limit the risks of oral floor lesion.

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Landes CA, Stübinger S, Ballon A, Sader R.

Piezoosteotomy in orthognathic surgery versus conventional saw and chisel osteotomy.

ENGLISH Oral Maxillofac Surg. 2008 Sep;12(3):139-47.

■ Piezoosteotomy in orthognathic surgery versus conventional saw and chisel osteotomy.

Introduction: Piezoosteotomy was assessed as alternative osteotomy method in orthognathic surgery regarding handling, time requirement, nerve and vessel impairment.

Materials and Methods: In this comparative clinical experience, 90 patient's orthognathic surgery procedures were performed in typical distribution prospectively by piezoosteotomy: 34 (38%) monosegment, 47 (52%) segmented LeFort I osteotomies, 94 (52%) sagittal split osteotomies, 11 (12%) symphyseal, and 4 (2%) mandibular body osteotomies. As controls served 90 retrospective patients with conventional saw and chisel osteotomy: 58 (64%) monosegment, 27 (30%) segmented LeFort I

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osteotomies, 130 (72%) sagittal split, and 4 (4%) symphyseal osteotomies.

Results and discussion: Piezoosteotomies were individually designed to interdigitate the jaw segments after repositioning. The pterygomaxillary suture weakened angulated tools; auxiliary chisels were required in 100% of cases for the nasal septum and lateral nasal walls, in 33% for pterygoid processes. The dorsal maxilla as the pterygoid process were easily reduced; 15% mandibular osteotomies required sawing, while the lingual dorsal osteotomy was performed by manual feedback due to limited visibility. Bloodloss decreased from average 537 +/- 208 ml vs. 772 +/- 338 ml ($p = 0.0001$). Operation time remained unchanged: 223 +/- 70 min vs. 238 +/- 60 min ($p = 0.2$) for a conventional bimaxillary procedure. Clinical courses and reossification were unobtrusive. Alveolar inferior nerve sensitivity was retained in 98% of the piezoosteotomy collective versus 84% of controls ($p = 0.0001$) at 3 months postoperative testing.

Conclusion: Piezoelectric osteotomy did not prolong the operation and reduced blood loss as alveolar nerve impairment. A few patients required additional sawing or chisel. Piezoelectric screw insertion as complex osteotomies may be initiated to simplify the procedure and increase segment interdigitation after repositioning as to minimize the osteofixation time and dimensions.

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■ **Intermaxillary fixation: technique and benefit for piezosurgical sagittal split osteotomy**

Introduction: The aim of this study was to assess piezosurgical sagittal split osteotomy with peroperative inter maxillary fixation.

Material and method: We studied 25 bimaxillary osteotomies, 50 sagittal split osteotomies performed with this technique. It included both maxillomandibular fixation during all the split osteotomy and performing split osteotomy in five steps. For each case, we noted the type of dysmorphia, the size of split osteotomy and the time required for surgery, along with common data such as sex, age, etc. The data was compared to results of a previous series of patients also operated with PIEZOSURGERY® but without peroperative maxillomandibular fixation.

Results: Using peroperative maxillomandibular fixation during piezosurgical bilateral sagittal osteotomy decreases the length of surgery by 33%, allows 9 times out of 10 for complete splitting, including the basilar edge, has no adverse effect especially on orthodontic material.

Discussion: PIEZOSURGERY® is a great progress for orthognathic surgery because of its precision and ability to preserve soft tissues. But it requires modification of the usual technique for mechanical section. Using peroperative inter maxillary fixation during ultrasonic splitting is a remarkably effective and easy technical modification.

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